



### **Product description**

Order number 65.716.321  
Product description Screw cap, blue, sterile, suitable for screw cap micro tubes, 500 piece(s)/double bag

### **Product characteristics**

Cap flat cap  
Closure type screw cap  
Suitable for screw cap micro tubes

### **Material & colours**

Product material Polypropylene (PP)  
Colour of product blue  
O-ring material Ethylene propylene diene rubber (EPDM)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Type of smallest subpackaging	double bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	540000
Depth of case	351 mm
Width of case	280 mm
Height of case	138 mm
Case volume	0.0136 cbm
Weight of case	2.05 kg
EAN of inner box	4038917923354
EAN case	4038917923361

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

