



**Product Specification**  
**Screw cap, pink, sterile, suitable for screw cap micro tubes**

**Page 1**



**Product description**

|                     |   |
|---------------------|---|
| Order number        | 65.716.324  |
| Product description | Screw cap, pink, sterile, suitable for screw cap micro tubes, 500 piece(s)/double bag |

**Product characteristics**

|              |                       |
|--------------|-----------------------|
| Cap          | flat cap              |
| Closure type | screw cap             |
| Suitable for | screw cap micro tubes |

**Material & colours**

|                   |  |
|-------------------|--|
| Product material  | Polypropylene (PP)                     |
| Colour of product | pink                                   |
| O-ring material   | Ethylene propylene diene rubber (EPDM) |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap, pink, sterile, suitable for screw cap  
micro tubes**

**Page 2**

**Purity & certification**

|                  |                                |
|------------------|--------------------------------|
| Product category | In vitro diagnostic device, CE |
| CE certified     | CE                             |
| Purity standard  | sterile                        |
| Sterilisation    | Electron irradiation           |
| Batched          | yes                            |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Type of smallest subpackaging | double bag    |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 5000          |
| Piece(s) / pallet             | 540000        |
| Depth of case                 | 351 mm        |
| Width of case                 | 280 mm        |
| Height of case                | 138 mm        |
| Case volume                   | 0.0136 cbm    |
| Weight of case                | 2.05 kg       |
| EAN of inner box              | 4038917923446 |
| EAN case                      | 4038917923453 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-02

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**