



Product description

Order number	65.716.328
Product description	Screw cap, violet, sterile, suitable for screw cap micro tubes, 500 piece(s)/double bag

Product characteristics

Cap	flat cap
Closure type	screw cap
Suitable for	screw cap micro tubes

Material & colours

Product material	Polypropylene (PP)
Colour of product	violet
O-ring material	Ethylene propylene diene rubber (EPDM)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Type of smallest subpackaging	double bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	540000
Depth of case	351 mm
Width of case	280 mm
Height of case	138 mm
Case volume	0.0136 cbm
Weight of case	2.05 kg
EAN of inner box	4038917923569
EAN case	4038917923576

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

