



Product description

| | |
|---------------------|---|
| Order number | 65.716.328 |
| Product description | Screw cap, violet, sterile, suitable for screw cap micro tubes, 500 piece(s)/double bag |

Product characteristics

| | |
|--------------|-----------------------|
| Cap | flat cap |
| Closure type | screw cap |
| Suitable for | screw cap micro tubes |

Material & colours

| | |
|-------------------|--|
| Product material | Polypropylene (PP) |
| Colour of product | violet |
| O-ring material | Ethylene propylene diene rubber (EPDM) |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap, violet, sterile, suitable for screw cap
micro tubes**

Purity & certification

| | |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD) |
| CE certified | CE - manufacturer's self-declaration |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Type of smallest subpackaging | double bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 5000 |
| Piece(s) / pallet | 540000 |
| Depth of case | 351 mm |
| Width of case | 280 mm |
| Height of case | 138 mm |
| Case volume | 0.0136 cbm |
| Weight of case | 2.05 kg |
| EAN of inner box | 4038917923569 |
| EAN case | 4038917923576 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

