



**Product Specification**  
**Push cap, natural, suitable for micro tube Ø 10.8 mm**

**Page 1**



**Product description**

Order number	65.723
Product description	Push cap, natural, suitable for micro tube Ø 10.8 mm, 1,000 piece(s)/bag

**Product characteristics**

Cap	standard cap
Closure type	push cap
Suitable for	micro tube Ø 10.8 mm

**Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

**Purity & certification**

Product category	In vitro diagnostic device, CE
CE certified	CE
Latex-free	yes
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-03-29

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Push cap, natural, suitable for micro tube Ø 10.8 mm**

**Page 2**

**Packaging**

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	640000
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0.0255 cbm
Weight of product	0.0004 kg
Weight of case	4.25 kg
EAN of inner box	4038917083713
EAN case	4038917002615

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-03-29

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**