



### **Product description**

Order number	72.689
Product description	Reaction tube, working volume: 2 ml, material: PP, transparent, cap: natural, retaining loop cap, cap attached, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Cap	retaining loop cap
Closure type	cap
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

### **Size**

Diameter	10.8 mm
Length of product	39 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	natural

### **Purity & certification**

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	150000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0014 kg
Weight of case	7.48 kg
EAN of inner box	4038917083416
EAN case	4038917002332

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

