



### **Product description**

Order number	72.690.378
Product description	Reaction tube, working volume: 1.5 ml, material: PP, yellow, cap: yellow, cap attached, with graduations, 500 piece(s)/bag

### **Product characteristics**

Graduation	yes
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g
Autoclavable (max. °C)	121 °C

### **Size**

Sample volume	1,5 ml
Length of product	39 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	yellow

### **Purity & certification**

Product category	no medical device   IVD
Batched	yes

### **Packaging**

Minimum order qty.	500
Type of inner packaging	bag
Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	160000
Case volume	0,0558 cbm
Weight of product	0,001 kg
Weight of case	5,6 kg
EAN of inner box	4038917369534
EAN case	4038917369527

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

