



### **Product description**

Order number	72.690.478
Product description	Reaction tube, working volume: 1.5 ml, material: PP, yellow, cap: yellow, retaining loop cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag

### **Product characteristics**

Graduation	yes
Cap	retaining loop cap
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g

### **Size**

Diameter	10.8 mm
Length of product	39 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow
Closure material	Polypropylene (PP)
Colour of cap	yellow

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	150000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.001 kg
Weight of case	5.69 kg
EAN of inner box	4038917161718
EAN case	4038917015257

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

