



Product description

Order number	72.695.200
Product description	SafeSeal reaction tube, working volume: 2 ml, material: PP, transparent, cap: transparent, SafeSeal cap, cap attached, with injected graduation and writing space, Biosphere® plus, 50 piece(s)/bag

Product characteristics

Graduation	yes
Cap	SafeSeal cap
Volume of work	2 ml
Centrifugation max (RCF)	25000 x g

Size

Diameter	10.8 mm
Length of product	40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	500
Piece(s) / pallet	72000
Depth of box	201 mm
Width of box	151 mm
Height of box	126 mm
Depth of case	323 mm
Width of case	213 mm
Height of case	145 mm
Case volume	0.01 cbm
Weight of product	0.0012 kg
Weight of case	0.94 kg
EAN of inner box	4038917117944
EAN case	4038917054256

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

