



**Product Specification**  
**SafeSeal reaction tube, 2 ml, PP, PCR**  
**Performance Tested**

**Page 1**



**Product description**

|                     |  |
|---------------------|--|
| Order number        | 72.695.400   |
| Product description | SafeSeal reaction tube, working volume: 2 ml, material: PP, transparent, cap: transparent, SafeSeal cap, cap attached, with injected graduation and writing space, PCR Performance Tested, 250 piece(s)/minigrip bag |

**Product characteristics**

|                          |              |
|--------------------------|--------------|
| Graduation               | yes          |
| Cap                      | SafeSeal cap |
| Volume of work           | 2 ml         |
| Centrifugation max (RCF) | 25000 x g    |

**Size**

|                   |         |
|-------------------|---------|
| Diameter          | 10.8 mm |
| Length of product | 40 mm   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**SafeSeal reaction tube, 2 ml, PP, PCR**  
**Performance Tested**

**Page 2**

**Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |
| Closure material  | Polypropylene (PP) |
| Colour of cap     | transparent        |

**Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Purity standard  | PCR Performance Tested        |
| Batched          | yes                           |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 2000          |
| Type of smallest subpackaging | minigrip bag  |
| Piece(s) / inner box          | 1000          |
| Piece (s) / outer case        | 2000          |
| Piece(s) / pallet             | 100000        |
| Depth of box                  | 301 mm        |
| Width of box                  | 221 mm        |
| Height of box                 | 161 mm        |
| Depth of case                 | 458 mm        |
| Width of case                 | 313 mm        |
| Height of case                | 180 mm        |
| Case volume                   | 0.0258 cbm    |
| Weight of product             | 0.0011 kg     |
| Weight of case                | 3.15 kg       |
| EAN of inner box              | 4038917223188 |
| EAN case                      | 4038917223171 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-02

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**