



### **Product description**

Order number	72.695.500
Product description	SafeSeal reaction tube, working volume: 2 ml, material: PP, transparent, cap: natural, SafeSeal cap, cap attached, with injected graduation and writing space, 250 piece(s)/bag

### **Product characteristics**

Graduation	yes
Cap	SafeSeal cap
Closure type	hinged lid
Volume of work	2 ml
Centrifugation max (RCF)	25000 x g

### **Size**

Diameter	10.8 mm
Length of product	40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	100000
Depth of box	301 mm
Width of box	221 mm
Height of box	161 mm
Depth of case	458 mm
Width of case	313 mm
Height of case	180 mm
Case volume	0.0258 cbm
Weight of product	0.0011 kg
Weight of case	3.09 kg
EAN of inner box	4038917266109
EAN case	4038917266093

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

