



### **Product description**

|                     |                                                                                                                        |
|---------------------|------------------------------------------------------------------------------------------------------------------------|
| Order number        | 72.699.002                                                                                                             |
| Product description | Reaction tube, working volume: 0.5 ml, material: PP, yellow, cap: yellow, standard cap, cap attached, 500 piece(s)/bag |

### **Product characteristics**

|                          |              |
|--------------------------|--------------|
| Graduation               | no           |
| Cap                      | standard cap |
| Closure type             | hinged lid   |
| Volume of work           | 0,5 ml       |
| Centrifugation max (RCF) | 20000 x g    |

### **Size**

|                   |        |
|-------------------|--------|
| Diameter          | 7.8 mm |
| Length of product | 30 mm  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | yellow             |
| Closure material  | Polypropylene (PP) |
| Colour of cap     | yellow             |

### **Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Batched          | yes                           |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 5000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 5000          |
| Piece(s) / pallet             | 320000        |
| Depth of case                 | 398 mm        |
| Width of case                 | 298 mm        |
| Height of case                | 215 mm        |
| Case volume                   | 0.0255 cbm    |
| Weight of product             | 0.0005 kg     |
| Weight of case                | 2.87 kg       |
| EAN of inner box              | 4038917118149 |
| EAN case                      | 4038917015899 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

