



Product description

Order number	72.699.003
Product description	Reaction tube, working volume: 0.5 ml, material: PP, blue, cap: blue, standard cap, cap attached, 500 piece(s)/bag

Product characteristics

Graduation	no
Cap	standard cap
Closure type	hinged lid
Volume of work	0,5 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	7.8 mm
Length of product	30 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	blue
Closure material	Polypropylene (PP)
Colour of cap	blue

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	320000
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0.0255 cbm
Weight of product	0.0005 kg
Weight of case	2.87 kg
EAN of inner box	4038917118156
EAN case	4038917015905

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

