



Product description

Order number	72.704.004
Product description	SafeSeal reaction tube, working volume: 0.5 ml, material: PP, yellow, cap: yellow, SafeSeal cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag

Product characteristics

Graduation	yes
Cap	SafeSeal cap
Volume of work	0,5 ml
Centrifugation max (RCF)	30000 x g

Size

Diameter	7.8 mm
Length of product	30 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	yellow
Closure material	Polypropylene (PP)
Colour of cap	yellow

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	156000
Depth of box	189 mm
Width of box	151 mm
Height of box	97 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0005 kg
Weight of case	1.61 kg
EAN of inner box	4038917227582
EAN case	4038917227575

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

