



Product Specification
SafeSeal reaction tube, 0.5 ml, PP, Biosphere® plus

Page 1



Product description

Order number	72.704.200
Product description	SafeSeal reaction tube, working volume: 0.5 ml, material: PP, transparent, cap: transparent, SafeSeal cap, cap attached, with injected graduation and writing space, Biosphere® plus, 50 piece(s)/bag

Product characteristics

Graduation	yes
Cap	SafeSeal cap
Volume of work	0,5 ml
Centrifugation max (RCF)	30000 x g

Size

Diameter	7.8 mm
Length of product	30 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
SafeSeal reaction tube, 0.5 ml, PP, Biosphere®
plus

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	500
Piece(s) / pallet	82500
Depth of box	189 mm
Width of box	151 mm
Height of box	97 mm
Depth of case	338 mm
Width of case	228 mm
Height of case	120 mm
Case volume	0.0092 cbm
Weight of product	0.0005 kg
Weight of case	0.62 kg
EAN of inner box	4038917235129
EAN case	4038917235112

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT