



### **Product description**

Order number	72.706.200
Product description	SafeSeal reaction tube, working volume: 1.5 ml, material: PP, transparent, cap: transparent, SafeSeal cap, cap attached, with injected graduation and writing space, Biosphere® plus, 50 piece(s)/bag

### **Product characteristics**

Graduation	yes
Cap	SafeSeal cap
Volume of work	1,5 ml
Centrifugation max (RCF)	30000 x g

### **Size**

Diameter	10.8 mm
Length of product	39 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**SafeSeal reaction tube, 1.5 ml, PP, Biosphere®**  
**plus**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	500
Piece(s) / pallet	72000
Depth of box	201 mm
Width of box	151 mm
Height of box	126 mm
Depth of case	323 mm
Width of case	213 mm
Height of case	145 mm
Case volume	0.01 cbm
Weight of product	0.0011 kg
Weight of case	0.9 kg
EAN of inner box	4038917118279
EAN case	4038917077491

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

