



Product description

Order number	72.706.201
Product description	SafeSeal reaction tube, working volume: 1.5 ml, material: PP, transparent, cap: transparent, SafeSeal cap, cap attached, with injected graduation and writing space, Biosphere® plus, 1 piece(s)/blister

Product characteristics

Graduation	yes
Cap	SafeSeal cap
Volume of work	1,5 ml
Centrifugation max (RCF)	30000 x g

Size

Diameter	10.8 mm
Length of product	39 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**SafeSeal reaction tube, 1.5 ml, PP, Biosphere®
plus**

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	120
Type of smallest subpackaging	blister
Piece(s) / inner box	60
Piece (s) / outer case	120
Piece(s) / pallet	17280
Depth of box	201 mm
Width of box	151 mm
Height of box	126 mm
Depth of case	323 mm
Width of case	213 mm
Height of case	145 mm
Case volume	0.01 cbm
Weight of product	0.001 kg
Weight of case	0.51 kg
EAN of inner box	4038917118286
EAN case	4038917079679

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT