



Product description

Order number	72.706
Product description	SafeSeal reaction tube, working volume: 1.5 ml, material: PP, transparent, cap: natural, SafeSeal cap, cap attached, with injected graduation and writing space, 250 piece(s)/bag

Product characteristics

Graduation	yes
Cap	SafeSeal cap
Volume of work	1,5 ml
Centrifugation max (RCF)	30000 x g

Size

Diameter	10.8 mm
Length of product	39 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	100000
Depth of box	301 mm
Width of box	221 mm
Height of box	161 mm
Depth of case	458 mm
Width of case	313 mm
Height of case	180 mm
Case volume	0.0258 cbm
Weight of product	0.001 kg
Weight of case	2.96 kg
EAN of inner box	4038917083577
EAN case	4038917073189

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

