



Product description

Order number	72.609.003
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, brown, without cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Knurls	yes
Base shape	skirted conical base
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product brown

Purity & certification

Satisfies the requirement IATA
Product category In-vitro diagnostic (IVD)
CE certified CE - manufacturer's self-declaration
Batched yes

Packaging

Minimum order qty. 5000
Type of smallest subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 5000
Piece(s) / pallet 180000
Depth of case 398 mm
Width of case 398 mm
Height of case 310 mm
Case volume 0.0491 cbm
Weight of product 0.0012 kg
Weight of case 6.86 kg
EAN of inner box 4038917114363
EAN case 4038917014267

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

