



Product description

| | |
|---------------------|---|
| Order number | 72.609.300 |
| Product description | Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, without cap, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------------|----------------------|
| Graduation | no |
| Knurls | yes |
| Base shape | skirted conical base |
| Volume of work | 2 ml |
| Centrifugation max (RCF) | 20000 x g |

Size

| | |
|-------------------|---------|
| Diameter | 10.8 mm |
| Length of product | 44 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Satisfies the requirement IATA
Product category In-vitro diagnostic, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 500
Type of smallest subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 5000
Piece(s) / pallet 180000
Case volume 0.0492 cbm
Weight of product 0.0012 kg
Weight of case 6.83 kg
EAN of inner box 4038917434287
EAN case 4038917434270

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

