



### **Product description**

Order number	72.609.601
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, without cap, with printed writing space, with graduations, 250 piece(s)/bag

### **Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Knurls	yes
Base shape	skirted conical base
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

### **Size**

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	IATA
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	5000
Piece(s) / pallet	180000
Depth of case	450 mm
Width of case	450 mm
Height of case	246 mm
Case volume	0.0498 cbm
Weight of product	0.0014 kg
Weight of case	7.06 kg
EAN of inner box	4038917795357
EAN case	4038917795364

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

