



Product description

Order number	72.693
Product description	Screw cap micro tube, working volume: 2 ml, conical base, with knurling, transparent, cap: natural, cap enclosed, 500 piece(s)/bag

Product characteristics

Graduation	no
Cap	standard screw cap
Closure type	screw cap
Knurls	yes
Base shape	conical base
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural
O-ring material	Ethylene propylene diene rubber (EPDM)

Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	5000
Type of inner packaging	bag
Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	180000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0015 kg
Weight of case	8.17 kg
EAN of inner box	4038917083454
EAN case	4038917002370

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

