



Product description

Order number	72.694.300
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, cap: natural, cap enclosed, 500 piece(s)/bag

Product characteristics

Graduation	no
Knurls	yes
Base shape	skirted conical base
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural
O-ring material	Ethylene propylene diene rubber (EPDM)

Purity & certification

Satisfies the requirement	IATA
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	160000
Depth of case	463 mm
Width of case	412 mm
Height of case	292 mm
Case volume	0.0557 cbm
Weight of case	8.59 kg
EAN of inner box	4038917926256
EAN case	4038917926263

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

