



### **Product description**

Order number	72.694.307
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, cap: natural, cap enclosed, with printed writing space, with graduations, 250 piece(s)/bag

### **Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	standard screw cap
Knurls	yes
Base shape	skirted conical base
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Diameter	10.8 mm
Length of product	44 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural
O-ring material	Ethylene propylene diene rubber (EPDM)

### **Purity & certification**

Satisfies the requirement	IATA
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	5000
Piece(s) / pallet	160000
Depth of case	463 mm
Width of case	412 mm
Height of case	292 mm
Case volume	0.0557 cbm
Weight of case	6.84 kg
EAN of inner box	4038917926386
EAN case	4038917926393

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

