



Product description

Order number	72.694.310
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, brown, cap: brown, cap enclosed, with printed writing space, with graduations, 250 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Knurls	yes
Base shape	skirted conical base
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Diameter	10.8 mm
Length of product	44 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	brown
Closure material	Polypropylene (PP)
Colour of cap	brown
O-ring material	Ethylene propylene diene rubber (EPDM)

Purity & certification

Satisfies the requirement	IATA
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	5000
Piece(s) / pallet	160000
Depth of case	463 mm
Width of case	412 mm
Height of case	292 mm
Case volume	0.0557 cbm
Weight of case	8.59 kg
EAN of inner box	4038917926423
EAN case	4038917926430

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

