



Product description

| | |
|---------------------|---|
| Order number | 72.694.700 |
| Product description | Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, cap: natural, cap assembled, with printed writing space, with graduations, PCR Performance Tested, Low DNA-binding, 100 piece(s)/minigrip bag |

Product characteristics

| | |
|--------------------------|----------------------|
| Label/ Print | with print |
| Colour of print/label | white |
| Graduation | yes |
| Cap | standard screw cap |
| Knurls | yes |
| Base shape | skirted conical base |
| Volume of work | 2 ml |
| Centrifugation max (RCF) | 20000 x g |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap micro tube, 2 ml, PCR Performance
Tested, Low DNA-binding

Page 2

Size

| | |
|-------------------|---------|
| Diameter | 10.8 mm |
| Length of product | 44 mm |

Material & colours

| | |
|-------------------|--|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | Polypropylene (PP) |
| Colour of cap | natural |
| O-ring material | Ethylene propylene diene rubber (EPDM) |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | IATA |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | minigrip bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 110000 |
| Depth of box | 211 mm |
| Width of box | 191 mm |
| Height of box | 141 mm |
| Depth of case | 398 mm |
| Width of case | 223 mm |
| Height of case | 160 mm |
| Case volume | 0.0142 cbm |
| Weight of product | 0.0016 kg |
| Weight of case | 2.11 kg |
| EAN of inner box | 4038917393447 |
| EAN case | 4038917393430 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-25

This document was prepared by EDP support and is valid without signature.



SARSTEDT