



Product Specification
Screw cap micro tube, 2 ml, sterile

Page 1



Product description

Order number	72.694.996
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, cap: colour mix, cap assembled, with printed writing space, with graduations, sterile, 100 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	standard screw cap
Knurls	yes
Base shape	skirted conical base
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Diameter	10.8 mm
Length of product	44 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	colour mix
O-ring material	Ethylene propylene diene rubber (EPDM)

Purity & certification

Satisfies the requirement	IATA
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	96000
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0.0135 cbm
Weight of product	0.0015 kg
Weight of case	1.86 kg
EAN of inner box	4038917117869
EAN case	4038917015721

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

