



Product description

Order number	72.703.004
Product description	Screw cap micro tubes, working volume: 1.5 ml, skirted conical base, with knurling, brown, without cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Knurls	yes
Base shape	skirted conical base
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	10,8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material Polypropylene (PP)
Colour of product brown

Purity & certification

Satisfies the requirement IATA
Product category In-vitro diagnostic, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 5000
Type of inner packaging bag
Piece(s) / inner box 500
Piece(s) / case 5000
Piece(s) / pallet 180000
Depth of case 398 mm
Width of case 398 mm
Height of case 300 mm
Case volume 0,0475 cbm
Weight of product 0,0013 kg
Weight of case 7,17 kg
EAN of inner box 4038917118217
EAN case 4038917071918

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

