



### **Product description**

Order number	72.703.416
Product description	Screw cap micro tube, working volume: 1.5 ml, skirted conical base, with knurling, transparent, cap: natural, cap attached assembled, with printed writing space, with graduations, PCR Performance Tested, 100 piece(s)/bag

### **Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	retaining loop cap
Knurls	yes
Base shape	skirted conical base
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap micro tube, 1.5 ml, PCR Performance Tested**

**Page 2**

**Size**

Diameter 10.8 mm  
Length of product 44 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent  
Closure material Polypropylene (PP)  
Colour of cap natural  
O-ring material Ethylene propylene diene rubber (EPDM)

**Purity & certification**

Satisfies the requirement IATA  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard PCR Performance Tested  
Batched yes

**Packaging**

Minimum order qty. 1000  
Type of smallest subpackaging bag  
Piece(s) / inner box 500  
Piece (s) / outer case 1000  
Piece(s) / pallet 110000  
Depth of box 211 mm  
Width of box 191 mm  
Height of box 141 mm  
Depth of case 398 mm  
Width of case 223 mm  
Height of case 160 mm  
Case volume 0.0142 cbm  
Weight of product 0.0017 kg  
Weight of case 2.27 kg  
EAN of inner box 4038917240659  
EAN case 4038917240642

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**