



**Product description**

Order number	72.703.700
Product description	Screw cap micro tube, working volume: 1.5 ml, skirted conical base, with knurling, transparent, cap: natural, cap assembled, with printed writing space, with graduations, PCR Performance Tested, Low DNA-binding, 100 piece(s)/minigrip bag

**Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	standard screw cap
Knurls	yes
Base shape	skirted conical base
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap micro tube, 1.5 ml, PCR Performance**  
**Tested, Low DNA-binding**

**Page 2**

**Size**

Diameter	10.8 mm
Length of product	44 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural
O-ring material	Ethylene propylene diene rubber (EPDM)

**Purity & certification**

Satisfies the requirement	IATA
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	110000
Depth of box	211 mm
Width of box	191 mm
Height of box	141 mm
Depth of case	398 mm
Width of case	223 mm
Height of case	160 mm
Case volume	0.0142 cbm
Weight of product	0.0016 kg
Weight of case	2.17 kg
EAN of inner box	4038917393232
EAN case	4038917393225

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-24

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**