



### **Product description**

|                     |   |
|---------------------|---|
| Order number        | 72.730.003  |
| Product description | Screw cap micro tube, working volume: 0.5 ml, skirted conical base, with knurling, transparent, without cap, 500 piece(s)/bag |

### **Product characteristics**

|                          |                      |
|--------------------------|----------------------|
| Graduation               | no                   |
| Knurls                   | yes                  |
| Base shape               | skirted conical base |
| Volume of work           | 0,5 ml               |
| Centrifugation max (RCF) | 20000 x g            |

### **Size**

|                   |         |
|-------------------|---------|
| Diameter          | 10.8 mm |
| Length of product | 44 mm   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |

### **Purity & certification**

|                           |                               |
|---------------------------|-------------------------------|
| Satisfies the requirement | IATA                          |
| Product category          | In vitrodiagnostic device, CE |
| CE certified              | CE                            |
| Batched                   | yes                           |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 5000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 5000          |
| Piece(s) / pallet             | 180000        |
| Depth of case                 | 398 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 300 mm        |
| Case volume                   | 0.0475 cbm    |
| Weight of product             | 0.0012 kg     |
| Weight of case                | 6.52 kg       |
| EAN of inner box              | 4038917119252 |
| EAN case                      | 4038917016520 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

