



Product description

Order number	72.730.100
Product description	Screw cap micro tube, working volume: 0.5 ml, skirted conical base, with knurling, transparent, cap: natural, cap attached, 500 piece(s)/bag

Product characteristics

Graduation	no
Cap	retaining loop cap
Knurls	yes
Base shape	skirted conical base
Volume of work	0,5 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural
O-ring material	Ethylene propylene diene rubber (EPDM)

Purity & certification

Satisfies the requirement	IATA
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	120000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0016 kg
Weight of case	8.7 kg
EAN of inner box	4038917119450
EAN case	4038917016643

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

