



### **Product description**

Order number	72.730.308
Product description	Screw cap micro tube, working volume: 0.5 ml, skirted conical base, with knurling, brown, cap: brown, cap enclosed, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Cap	standard screw cap
Closure type	screw cap
Knurls	yes
Base shape	skirted conical base
Volume of work	0,5 ml
Centrifugation max (RCF)	20000 x g

### **Size**

Diameter	10.8 mm
Length of product	44 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	brown
Closure material	Polypropylene (PP)
Colour of cap	brown
O-ring material	Ethylene propylene diene rubber (EPDM)

### **Purity & certification**

Satisfies the requirement	IATA
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

### **Packaging**

Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	150000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap micro tube, 0.5 ml**

**Page 2**

Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of case	8.23 kg
EAN of inner box	4038917994408
EAN case	4038917994415

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2026-03-08

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**