



Product description

Order number	72.730.309
Product description	Screw cap micro tube, working volume: 0.5 ml, skirted conical base, with knurling, transparent, without cap, with printed writing space, 250 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	no
Knurls	yes
Base shape	skirted conical base
Volume of work	0,5 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	44 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	5000
Piece(s) / pallet	150000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap micro tube, 0.5 ml

Page 2

Case volume	0.0475 cbm
Weight of case	6.62 kg
EAN of inner box	4038917994484
EAN case	4038917994491

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-04-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT