



### **Product description**

Order number	72.785
Product description	Screw cap micro tube, working volume: 0.5 ml, skirted conical base, transparent, without cap, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Knurls	no
Base shape	skirted conical base
Volume of work	0,5 ml
Centrifugation max (RCF)	20000 x g

### **Size**

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	IATA
Product category	In-vitro diagnostic, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	5000
Type of inner packaging	bag
Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	180000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0012 kg
Weight of case	6.63 kg
EAN of inner box	4038917084116
EAN case	4038917077675

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

