



Product description

Order number	79.8701.502
Product description	Blood bag, nominal volume: 500 ml, with venous cannula, preparation: CPDA-1 stabilizer, with cannula protection, sterile, 1 piece(s)/bag

Product characteristics

Type of preparation	CPDA-1 stabilizer
Application	Obtaining blood and blood products

Size

Sample volume	500 ml
Width of product	21 cm
Length of product	27 cm

Material & colours

Colour of product	transparent
-------------------	-------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category	Medical device
CE certified	CE 0123 (TÜV SÜD)
Purity standard	sterile
Pyrogen-free	yes
Plasticiser	Diethylhexylphthalate (DEHP)
Batched	yes

Packaging

Minimum order qty.	40
Type of inner packaging	bag
Piece(s) / inner box	40
Piece(s) / case	40
Case volume	0.0311 cbm
Weight of product	0.1258 kg
Weight of case	5.6 kg
EAN of inner box	4038917148450
EAN case	4038917031424

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

