



Product description

Order number	79.8701.990
Product description	Cap, for blood bags, to seal the neck on blood bags

Product characteristics

Application	Obtaining blood and blood products
Suitable for	blood bags

Material & colours

Colour of product	red
-------------------	-----

Purity & certification

Product in accordance with MDD or IVDD	Medical device
CE certified	CE0482 (Medcert)
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Lot or serial number	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	700
Piece(s) / inner box	700
Piece(s) / case	700
Piece(s) / pallet	44800
Case volume	0,0199 cbm
Weight of product	0,0017 kg
Weight of case	1,57 kg
EAN of inner box	4028339690224
EAN case	4028339690217

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

