



Product description

Order number	88.7750.415
Product description	RadTag® irradiation indicator X-ray, irradiation factor: 15/50 Gy, for visual monitoring of the minimum and maximum irradiation intensity of blood and blood products

Product characteristics

Application	Transfusion and infusion
-------------	--------------------------

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	200
Type of smallest subpackaging	case
Piece(s) / inner box	200
Piece (s) / outer case	200
Piece(s) / pallet	64800
Depth of case	268 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0004 kg
Weight of case	0.22 kg
EAN of inner box	4038917260374
EAN case	4038917260367

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

