



Product description

Order number	88.7780.190
Product description	Box container, (LxWxH): 95 x 60 x 60 mm, for storing blood products

Product characteristics

Application	Transfusion and infusion
-------------	--------------------------

Size

Width of product	60 mm
Height of product	60 mm
Length of product	95 mm

Material & colours

Colour of product	white
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 300
Piece(s) / inner box 1
Piece(s) / case 300
Piece(s) / pallet 14400
Case volume 0,0245 cbm
Weight of product 0,011 kg
Weight of case 3,72 kg
EAN of inner box 4038917142014
EAN case 4038917080767

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

