



Product description

Order number	88.7780.197
Product description	Box small, (LxWxH): 200 x 130 x 25 mm, with window, for storing blood products

Product characteristics

Application	Transfusion and infusion
Viewing window	yes

Size

Width of product	130 mm
Height of product	25 mm
Length of product	200 mm

Material & colours

Colour of product	white
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 300
Type of inner packaging case
Piece(s) / inner box 1
Piece(s) / case 300
Piece(s) / pallet 8400
Case volume 0.0502 cbm
Weight of product 0.029 kg
Weight of case 9.42 kg
EAN of inner box 4038917142069
EAN case 4038917063159

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

