



Product description

Order number 95.1008
Product description TourniSTRETCH tourniquet, latex-free, 25 piece(s)/bag

Product characteristics

Application venous blood collection

Material & colours

Colour of product yellow

Purity & certification

Product in accordance with MDD or IVDD Medical device
CE certified CE supplier > no Sarstedt CE
Latex-free yes
Lot or serial number yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

| | |
|-------------------------|---------------|
| Minimum order qty. | 25 |
| Type of inner packaging | bag |
| Piece(s) / inner box | 25 |
| Piece(s) / case | 500 |
| Piece(s) / pallet | 8000 |
| Case volume | 0,0572 cbm |
| Weight of product | 0,0087 kg |
| Weight of case | 5,2 kg |
| EAN of inner box | 4038917405676 |
| EAN case | 4038917405669 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

