



Product description

Order number 95.1008
Product description TourniSTRETCH tourniquet, latex-free, 25 piece(s)/bag

Product characteristics

Application venous blood collection

Material & colours

Colour of product yellow

Purity & certification

Product in accordance with MDD or IVDD Medical device
CE certified CE supplier > no Sarstedt CE
Latex-free yes
Lot or serial number yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	25
Type of inner packaging	bag
Piece(s) / inner box	25
Piece(s) / case	500
Piece(s) / pallet	8000
Case volume	0,0572 cbm
Weight of product	0,0087 kg
Weight of case	5,2 kg
EAN of inner box	4038917405676
EAN case	4038917405669

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

