



Product description

Order number	97.8700.330
Product description	Mobile blood donation table, Side table for mobile blood collection, can be collapsed quickly and easily, (LxWxH): 40 x 55 x 73 cm

Product characteristics

Application	Obtaining blood and blood products
-------------	------------------------------------

Size

Width of product	55 cm
Height of product	73 cm
Length of product	40 cm

Material & colours

Colour of product	red
-------------------	-----

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product in accordance with MDD or IVDD	no medical device IVD
Lot or serial number	no

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Piece(s) / pallet	10
Case volume	0,1014 cbm
Weight of product	6 kg
Weight of case	7 kg
EAN of inner box	4028339671421
EAN case	4028339671414

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

