



Product description

Order number	97.8700.452
Product description	Software, for backing up donation and registration data on a host computer and for central donor monitoring

Product characteristics

Application	Obtaining blood and blood products
-------------	------------------------------------

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Weight of product	0,0005 kg
EAN of inner box	4028339633924
EAN case	4028339633917

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

