



## Product description

Order number	97.8700.980
Product description	Precision weight, for weight calibration

## Product characteristics

Application	Obtaining blood and blood products
-------------	------------------------------------

## Material & colours

Colour of product	blue
-------------------	------

## Purity & certification

Product in accordance with MDD or IVDD	no medical device   IVD
Lot or serial number	no

## Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0,0003 cbm
Weight of product	0,5 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification  
Calibration weight 500 g**

**Page 2**

Weight of case	0,54 kg
EAN of inner box	4028339634822
EAN case	4028339634815

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-07-24

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**