



Product description

Order number	97.8700.423
Product description	TS-Power-Pack Li-Ion, Li-ion battery for mains-independent operation of the TOPSWING, with integrated charging unit

Product characteristics

Application	Obtaining blood and blood products
Battery power	6750 mAh

Size

Width of product	7 cm
Height of product	7 cm
Length of product	20 cm

Material & colours

Colour of product	black
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category	Medical device
CE certified	CE
Batched	no

Packaging

Minimum order qty.	1
Type of inner packaging	case
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0.0039 cbm
Weight of product	0.95 kg
Weight of case	1.09 kg
EAN of inner box	4028339632422
EAN case	4028339632415

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

