



### **Product description**

Order number	97.8700.432
Product description	Module TS-Transport, Storage and transport cover for TOPSWING with height adjustment by approx. 9 cm during donation, stackable

### **Product characteristics**

Application	Obtaining blood and blood products
-------------	------------------------------------

### **Size**

Width of product	26 cm
Height of product	10 cm
Length of product	48 cm

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





## **Packaging**

Minimum order qty.	1
Type of inner packaging	case
Piece(s) / inner box	1
Piece(s) / case	1
Piece(s) / pallet	56
Case volume	0.03 cbm
Weight of product	1 kg
Weight of case	1.5 kg
EAN of inner box	4028339632927
EAN case	4028339632910

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

