



Product description

Order number	97.8700.433
Product description	Module TS-Level, for adjusting the height of the TOPSWING by approx. 9 cm during donation, stackable

Product characteristics

Application	Obtaining blood and blood products
-------------	------------------------------------

Size

Width of product	26 cm
Height of product	10 cm
Length of product	48 cm

Purity & certification

Product category	no medical device IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0,03 cbm
Weight of product	1 kg
Weight of case	1,5 kg
EAN of inner box	4028339633023
EAN case	4028339633016

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

