



Product description

Order number 97.8700.433
Product description Module TS-Level, for adjusting the height of the TOPSWING by approx. 9 cm during donation, stackable

Product characteristics

Application Obtaining blood and blood products

Size

Width of product 26 cm
Height of product 10 cm
Length of product 48 cm

Purity & certification

Product category no medical device | IVD
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	1
Type of inner packaging	case
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0.03 cbm
Weight of product	1 kg
Weight of case	1.5 kg
EAN of inner box	4028339633023
EAN case	4028339633016

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

