



Product description

Order number	97.8700.440
Product description	Scanner TOUCH 65 PRO, Barcode scanner with connection for the TOPSWING

Product characteristics

Application	Obtaining blood and blood products
-------------	------------------------------------

Size

Width of product	7 cm
Length of product	16 cm

Material & colours

Colour of product	grey
-------------------	------

Purity & certification

Product category	Appliance, according to machinery or equipment directive, CE
CE certified	CE
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	1
Type of inner packaging	case
Piece(s) / inner box	1
Piece(s) / case	1
Piece(s) / pallet	260
Case volume	0.0056 cbm
Weight of product	0.24 kg
Weight of case	0.39 kg
EAN of inner box	4028339633429
EAN case	4028339633412

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

